220.00 mm

158.00 mm

RENNOVA[®]

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COMPOSITION:
Highly cross-linked sodium hyaluronate at 23mg/mL with 0.32% lidocaine hydrochloride. Phosphate buffer pH 7.3 q.s.

DESCRIPTION:
RENNOVA LIDO[®] is a sterile, transparent viscoelastic gel containing 23mg/mL of sodium hyaluronate. Sodium hyaluronate (sodium salt of hyaluronic acid) belongs to a group of few identical substances in all living organisms. It is the natural polysaccharide present in all body tissues, which have particularly high concentrations in synovial fluid and skin. It is a single-use product.
Model: Rennova Lift Ultra Volume Lido 1mL | Rennova Lift Ultra Volume Lido 2mL.

INDICATIONS:
The **RENNOVA LIDO[®]** device is intended for the treatment of deep wrinkles on the face resulting from the intrinsic and extrinsic aging process and is therefore recommended for use on puppet lines and nasolabial folds. The maximum recommended dosage is 10mL per application and 20mL per year, in medium to deep dermis. However, the maximum volumes to be applied in each of the areas are defined by qualified health professionals according to the patient's previous evaluation and the condition of the skin texture and/or severity of wrinkles and expression lines. It is recommended to use it in the following areas of application, leaving it up to the physician to choose which one is the most suitable for the patient. Recommended application areas: Puppet Lines and Nasolabial Groove.

CONTRAINDICATIONS
1- For patients with a history of anaphylaxis or severe allergic symptoms.
2- For pregnant or breastfeeding women.
3- The use is not recommended in areas where other products are already being used.
4- For patients who are or were hypersensitive to hyaluronic acid.
5 - For patients who are or have been hypersensitive to amide-type local anesthetics such as lidocaine.
6- For patients who are or were hypersensitive to the raw materials of the filling material.
7- For patients with bleeding disorder and any autoimmune disease.
8- For correction of superficial wrinkles.
9- For combination with peeling or laser treatment.
10 - For patients with skin disorder and inflammation at the area under treatment.
11 - Do not use the product in blood vessels, use the product carefully.
Contraindicated Areas: Glabellar or periorcular area; Periorbital area; Perioral region; superficial dermis; Injection into vessels; Injection into the epidermis.

PRECAUTIONS
Normal precautions associated with intradermal injections must be taken. Like any procedure of this type, the **RENNOVA LIDO[®]** implant is associated with an inherent risk of infection. **RENNOVA LIDO[®]** must be used on or near anatomical areas where there is active skin disease, inflammation, or related conditions. In patients who are taking aspirin or non-steroidal anti-inflammatory drugs, injection of the product may increase the severity of bleedings at the injection area and must be carefully injected in unavoidable situations. Do not inject it intravascularly. **RENNOVA LIDO[®]** has not been tested on pregnant women or children. The product must be checked before use and must not be used if its conditions are not the regular ones.

"SINGLE-USE PRODUCT. REPROCESSING PROHIBITED"

CAUTIONS
The qualified healthcare professional must carefully follow the general precautions regarding the hypodermic area and the superior periosteum one. It is advisable to aspirate before injecting the implant and to take special care in highly vascularized areas. The product must not be injected into the blood vessels, which may result in occlusion, ischemia, infarct, and necrosis of local and distant tissues.
1 - Operating surgeons must have sufficient anatomical knowledge about the area in which the surgery is being carried out.
2 - Operating surgeons must always be careful about the risk of contamination of operating rooms and equipment used for the procedure.
3 - Operating surgeons must carefully follow general precautions regarding the hypodermic area and the superior periosteum one.
4 - Operating surgeons must not inject the product into areas with permanent implants.
5 - Operating surgeons must examine whether the products remain sterile before use or not.
6 - Operating surgeons must check the expiration date on the product label.
7 - All remaining content of the products must not be reused and must be discarded.
8 - If the needle package is by chance broken, the needle must not be used. Injection needles must be discarded.
9 - The product must not be used for patients who anticipate more than specified.
10 - Patients must not put on makeup for 12 hours after the injection. The injection area shall not be exposed to sunlight, UV rays and high and low temperature environments such as a sauna for 2 weeks after the injection.
11 - If pain is detected while the product is being injected, the injection must be stopped immediately. If pain is detected after the injection, consult a specialist about this symptom.
12 - Operating surgeons must perform the treatment after fully understanding the instructions for use.

INCOMPATIBILITY
As hyaluronic acid can be precipitated by quaternary ammonium compounds and chlorhexidine, such as specific chemical sterilizers, it must not be stored near these materials or handled with operating equipment that came into contact with these materials.

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220.00 mm

158.00 mm

Redness, swelling, sensitiveness, pain, bruising, itching, and other symptoms may occur in the surgical area, which will naturally disappear in a week or two. In such cases, oral administration of corticosteroids may be effective. Patients who experience these symptoms cannot be treated with **RENNOVA LIDO**® or other similar products. Operating surgeons shall inform patients that the following symptoms occur immediately or later after injection and must report the occurrence of these symptoms.

- After the injection, symptoms of inflammation (fever, redness, swelling, etc.) with soreness or pain may persist for about a week.
- Be careful about the possibility of abscess or hypersensitivity occurring after injection.
- If the inflammation symptom persists for more than a week or another adverse event occurs, consult a specialist immediately.

HOW TO USE

The product must only be injected by a doctor or qualified health professional in accordance with current legislation in the country. Before injecting a doctor must examine the patient's physical condition. If the syringe blister pack is open or damaged, do not use it. Do not mix **RENNOVA LIDO**® with other injectable materials. **RENNOVA LIDO**® must be injected using a strict aseptic technique. Remove the cap from the syringe tip and attach an appropriate sterilized needle: the recommended size is 21G (gauge) to 30G – (The kit does not come with a needle or microcannula). The injection area must be cleaned with alcohol or other suitable antiseptic solution before injection. Inject slowly, the injection volume depends on the depth of the wrinkles. Massage the injected area. For the correct conditions to persist, additional periodic injections of the product may be required, and the doctor must inform the patient of that.

Note: The kit is provided without needle or microcannula.

PRESENTATION

RENNOVA LIDO® is supplied in a pre-filled plastic syringe, blister pack, in the following models:
1 syringe with 1mL of Rennova Lift Ultra Volume Lido;
1 syringe with 2mL of Rennova Lift Ultra Volume Lido.
The contents of the syringe are sterile and endotoxin-free. Store and transport it between 2°C and 25°C. Do not freeze it.
Conditions: Do not expose it directly to sunlight.

EXPIRATION DATE

The expiration date is indicated on the package.

DEGRADATION TIME

The average period of absorption/degradation of the product after implantation is 6 to 12 months. However, it can vary according to the metabolism of each individual.

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HUMEDIX CO., LTD.